

Next Generation Knee Instruments

**SAYAN**<sup>®</sup>  
ORTHOPAEDICS

Harmony

effectiveness



info@sayan.com.tr  
sayan.com.tr

**Sayan Orthopaedics**

Tel : +90 232 472 27 78 Fax : +90 232 472 27 76



## External Reference Tibial Cutting Guide

- 1 10° cutting slope angle can be modified to the design of the implant
- 2 Right, left and neutral versions are available for the tibial cutting guide
- 3 Spiked and non-spiked fixation versions are available according to Surgeon's preference
- 4 Length of the system is adjusted to patient's by locking cam
- 5 The system could be fixed around the ankle either by ankle clamps or by using spring or both for more rigid fixation when needed
- 6 Fine M/L positioning of the system by distal knob
- 7 Alignment of the guide for posterior sloped cut

