

Next Generation Knee Instruments

SAYAN[®]
ORTHOPAEDICS


Sturdiness



info@sayan.com.tr
sayan.com.tr

Sayan Orthopaedics

Tel : +90 232 472 27 78 Fax : +90 232 472 27 76



Patella Saw Guide

- 1** Minimum of 10 mm resection is adjusted according to design of the implant
- 2** Single-button resection adjustment
- 3** Graduated increments depending on the size of the implant
- 4** Saw slot can be adjusted to any thickness
- 5** Sharp spikes to hold the patella tight
- 6** Position retaining by the locking nut on the side
- 7** Saw slot with extended base for easier aiming of the blade

