

Customer-made
Precision Broaching



Broach

- Modified according to the coupling of the broach handle
- Precisely grinded to the section of the corresponding implant
- Variable pitch of cutting teeth

fine art




Next Generation Orthopaedic Instruments

SAYAN[®]
ORTHOPAEDICS




excellence



Commonly machined sharp edges :

- Less durability
- Smaller bone cavity as used
- Easy to get damaged therefore shorter life time



SAYAN double cutting angle :

- More durable with R0.4 mm on the cutting tip
- Better bone engraving with double cutting angle
- Longer life time



sayan@sayan.com.tr
sayan.com.tr

Sayan Orthopaedics

Tel : +90 232 472 27 78 Fax : +90 232 472 27 76