

Flexible Solutions for MIS Reamering



MIS-**Flex** Acetabular Reamer Handle

sayan@sayan.com.tr sayan.com.tr Sayan Orthopaedics

Tel: +90 232 472 27 78 Fax: +90 232 472 27 76



MIS-Flex Acetabular Reamer Handle

6

7
2

- Modified according to coupling of the reamer
- 2 Hinge to adjust the handle for left-hand and right-hand usage of the motor, sweeping 86° in total
- 3 Ergonomic handle to hold and drive the instrument
- Quick coupling end available in various motor coupling types

- 5 Provides 10° extra angular freedom to drive the instrument, enabling the handle sweep 106° in total from right to left
- Protective cover to avoid rotating shaft pining gloves
- 7 Alternative handle to hold and drive the instrument in case of free-hand use by disassembling the ergonomic handle at the back. During free-hand use, the motor driven shaft has a freedom of motion upto 81° vertically and 67° horizontally on one side



Sayan Orthopaedics

sayan@sayan.com.tr sayan.com.tr

Tel: +90 232 472 27 78

Fax: +90 232 472 27 76