

Next Generation Knee Instruments

SAYAN[®]
ORTHOPAEDICS

Sturdiness



sayan@sayan.com.tr
sayan.com.tr

Sayan Orthopaedics

Tel : +90 232 472 27 78 Fax : +90 232 472 27 76



Patella Saw Guide

- 1 Minimum of 10 mm resection is adjusted according to design of the implant
- 2 Single-button resection adjustment
- 3 Graduated increments depending on the size of the implant
- 4 Saw slot can be adjusted to any thickness
- 5 Sharp spikes to hold the patella tight
- 6 Position retaining by the locking nut on the side
- 7 Saw slot with extended base for easier aiming of the blade

