


# Next Generation Spinal Instruments

**SAYAN**<sup>®</sup>  
ORTHOPAEDICS

Easy hold  
handle



strength

info@sayan.com.tr  
sayan.com.tr

**Sayan Orthopaedics**

Tel : +90 232 472 27 78 Fax : +90 232 472 27 76

## Lateral Persuader

- 1 Tip to be modified according to design of the tulip
- 2 Self-locking of the tulip by screwing-in the persuading shaft
- 3 Perforated all around for efficient cleaning
- 4 Comfort grip handle
- 5 Position retaining handle



rigidity

